

| BRADENHEAD INSPECTION REPORT | | | |
|--|---------------------------------|---------------|-------------|
| Stop Name | WEAVER RIDGE #23-14 - CO085-009 | Inspected By: | Jacob White |
| Inspection Date | 9/10/24 | Timestamp | 3:59 PM |
| Bradenhead Test | | | |
| Well Status | | | Status |
| Number of Casing Strings | | | Flowing |
| Multiple Completions | | | Two |
| Existing Pressures | | | No |
| Formation Tubing (Existing) | | | 37 |
| Formation Tubing -2 (Existing) | | | |
| Production Casing | | | 71 |
| Intermediate Casing | | | |
| Surface Casing | | | 0 |
| Elapsed Time 0 Minutes | | | |
| Formation Tubing | | | 37 |
| Formation Tubing -2 | | | |
| Production Casing PSIG | | | 70 |
| Intermediate Casing PSIG | | | |
| Bradenhead Flow | | | No Flow |
| Bradenhead Fluid | | | Vapor |
| Bradenhead Fluid Description at 0 Minutes | | | -- |
| Elapsed Time 5 Minutes | | | |
| Formation Tubing - 5 Min | | | 37 |
| Formation Tubing -2 - 5 Min | | | |
| Production Casing PSIG - 5 Min | | | 71 |
| Intermediate Casing PSIG - 5 Min | | | |
| Bradenhead Flow - 5 Min | | | No Flow |
| Bradenhead Fluid - 5 Min | | | None |
| Bradenhead Fluid Description at 5 Minutes | | | -- |
| Elapsed Time 10 Minutes | | | |
| Formation Tubing - 10 Min | | | 37 |
| Formation Tubing -2 - 10 Min | | | |
| Production Casing PSIG - 10 Min | | | 71 |
| Intermediate Casing PSIG - 10 Min | | | |
| Bradenhead Flow - 10 Min | | | No Flow |
| Bradenhead Fluid - 10 Min | | | None |
| Bradenhead Fluid Description at 10 Minutes | | | -- |
| Elapsed Time 15 Minutes | | | |
| Formation Tubing - 15 Min | | | 37 |
| Formation Tubing -2 - 15 Min | | | |
| Production Casing PSIG - 15 Min | | | 72 |
| Intermediate Casing PSIG - 15 Min | | | |
| Bradenhead Flow - 15 Min | | | No Flow |
| Bradenhead Fluid - 15 Min | | | None |
| Bradenhead Fluid Description at 15 Minutes | | | -- |
| Elapsed Time 20 Minutes | | | |
| Formation Tubing - 20 Min | | | 37 |
| Formation Tubing -2 - 20 Min | | | |
| Production Casing PSIG - 20 Min | | | 72 |
| Intermediate Casing PSIG - 20 Min | | | |
| Bradenhead Flow - 20 Min | | | No Flow |
| Bradenhead Fluid - 20 Min | | | None |
| Bradenhead Fluid Description at 20 Minutes | | | -- |
| Elapsed Time 25 Minutes | | | |
| Formation Tubing - 25 Min | | | 36 |
| Formation Tubing -2 - 25 Min | | | |
| Production Casing PSIG - 25 Min | | | 69 |
| Intermediate Casing PSIG - 25 Min | | | |
| Bradenhead Flow - 25 Min | | | No Flow |
| Bradenhead Fluid - 25 Min | | | None |
| Bradenhead Fluid Description at 25 Minutes | | | -- |
| Elapsed Time 30 Minutes | | | |
| Formation Tubing - 30 Min | | | 36 |
| Formation Tubing -2 - 30 Min | | | |
| Production Casing PSIG - 30 Min | | | 73 |
| Intermediate Casing PSIG - 30 Min | | | |
| Bradenhead Flow - 30 Min | | | No Flow |
| Bradenhead Fluid - 30 Min | | | None |
| Bradenhead Fluid Description at 30 Minutes | | | -- |
| Elapsed Time 5 Minutes SI | | | |
| Formation Tubing - 5 Min SI | | | 36 |
| Formation Tubing -2 - 5 Min SI | | | |
| Production Casing PSIG - 5 Min SI | | | 74 |
| Intermediate Casing PSIG - 5 Min SI | | | |
| Bradenhead Flow - 5 Min SI | | | No Flow |
| Bradenhead Fluid - 5 Min SI | | | None |
| Bradenhead Fluid Description at 5 Minutes SI | | | -- |
| Instantaneous Bradenhead Pressure at End of Test | | | 0 |
| Buried Valve? | | | No |
| Confirmed Open? | | | Yes |
| Bradenhead Sample Taken? | | | No |
| Bradenhead Sample Type | | | -- |
| Note Size of Valve: BHD | | | 2" |
| Note Size of Valve: INTERMEDIATE | | | 1/2" |
| | | | |
| Remarks | null | | |